

Current information on Rapid Diagnostic Tests

Rapid Tests for HIV Commercially Available Rapid Tests for HIV

We encourage prospective buyers of rapid tests to request and review all available information and to make their own judgement in selecting tests most suitable for their specific programs. The following list of known suppliers is intended only to provide a starting point for these enquiries. Also, note that many of these manufacturers have distribution capabilities in a number of countries.

This list is also available for download in Rich Text Format.

- [Download](#) (70 Kb .rtf)

Time: All tests included in these tables can be completed in approximately 30 minutes or less.

Equipment: Tests may need blood draw material, centrifuge, pipettes, pipette tips, and specimen containers for serum- and plasma-based tests. For tests that use whole blood, test tubes, test tube racks, and capillary tubes may be needed.

Type of Test: Lateral Flow

Manufacturer	Name of Test	Web Site	Useable Specimen
Abbott Diagnostics	Determine HIV-1/2	http://www.abbottdiagnostics.com	Plasma, serum, whole blood
Akers Laboratories	Health Test HIV-1 Assay	http://www.akerslaboratories.com	
AmeriTek, Inc.	dBest one-step HIV-	http://www.ameritek.org	Serum, whole

	1/HIV-2		blood
Biokit, S.A.	Biorapid HIV 1&2	http://www.biokit.com	Serum, plasma, whole blood
BioRad Laboratories	Genie II HIV 1/2	http://www.biorad.com	Serum, plasma
ChemBio Diagnostics Systems, Inc.	HIV 1/2 Stat-Pak	http://www.chembio.com	Plasma, serum, whole blood
Cypress Diagnostics	HIV 1/2 Quick Test	http://www.diagnostics.be/cypress	Serum
International Newtech Development	One Step HIV (1+2) Serum Test	http://www.ind.ca	Serum
Makro Medical (Pty) Ltd.	HIV 1/2 rapid test	http://users.iafrica.com/m/ma/makromed/id3.htm	Whole blood
Mega Diagnostics, Inc.	MegaKwik HIV 1&2 Card Test	http://www.mega-dx.com	
Orgenics	ImmunoGold HIV &2	http://www.orgenics.com	Serum, plasma
OraSure Technologies, Inc.	OraQuick HIV-1/2	http://www.epitope.com	Saliva
Princeton BioMeditech Corporation	BioSign HIV	http://www.pbmc.com	
Genix Technology	RapidTech Plus HIV 1+2	http://www.genixtech.com	Serum, whole blood
Saliva Diagnostic System	Hema Strip HIV	http://www.salv.com	Whole Blood
Saliva Diagnostic System	Saliva Strip HIV	http://www.salv.com	Saliva
Saliva Diagnostic System	Sero Strip HIV	http://www.salv.com	Serum, plasma
Standard Diagnostics, Inc.	HIV- 1/2 3.0	http://www.standardia.com	Serum, plasma, whole blood

Trinity Biotech	Uni-Gold HIV Peptide	http://www.trinitybiotech.com	Serum, plasma, whole blood
Trinity Biotech	Uni-Gold HIV Recombinant	http://www.trinitybiotech.com	Serum, plasma, whole blood
World Diagnostics, Inc.	Smart Check HIV 1&2 Whole Blood	http://www.globalemed.com	Serum, plasma, whole blood
World Diagnostics, Inc.	Smart Check HIV 1&2 Double Spot	http://www.globalemed.com	Serum, plasma

Type of Test: Flow-Through

Manufacturer	Name of Test	Web Site	Use & Spec
BioRad Laboratories	Multispot	http://www.biorad.com	Serum Plasma
Diagnostics Product Source	Rapid HIV 1/2 flow through spot test	http://www.aptechservices.com/diag/infectious.html	Serum plasma
EY Laboratories, INC.	HIV SPOT-1/2 / HIV CHEK 1/2	http://www eylabs.com/	Serum plasma
J. Mitra & Co. Ltd.	HIV Tri-Dot	http://www.jmltd.com	Plasma serum
MedMira Laboratories	MedMira HIV-1/2	http://www.medmira.com	Serum plasma whole blood
Organics	DOUBLECHECK HIV1&2	http://www.orgenics.com	Serum plasma
Sayvon Diagnostics Ltd (Diatech)	HIV SAV 1/2	http://hctech.com/sayvon/pages/main.html	Serum plasma
Specialty Systems, Inc.	AccuSpot Plus HIV 1&2	http://www.accudx.com	Serum plasma
Trinity Biotech	SeroCard HIV	http://www.trinitybiotech.com	Serum plasma whole blood

Trinity Biotech	SalivaCard HIV	http://www.trinitybiotech.com	Saliv
UNIVERSAL Healthwatch	QUIXTM HIV 1 (M&O), HIV 2 RAPID BLOOD TEST DEVICE	http://www.guardianscientificafrika.com	Who Bloo

Type of Test: Agglutination

Manufacturer	Name of Test	Web Site	Useable Specimen	Perf. Ref.
Fujirebio	Serodia HIV-1/2	http://www.fujirebioamerica.com	Serum, plasma	(4), (5), (16), (20)
Trinity Biotech	Capillus HIV	http://www.trinitybiotech.com	Serum, plasma, whole blood	(4), (5), (9), (16), (18), (20)

Unknown

Manufacturer	Name of Test	Web Site	Useable Specimen	Perf. Ref.
AMDL, Inc.	Q Spot HIV 1&2 Components	http://www.amdl.com		

Type of Test: Magnetic Beads

Manufacturer	Name of Test	Web Site	Useable Specimen	Performance References
Bionor A/S	Bionor HIV-1/2	http://www.bionor.no	Serum, plasma, whole blood	(6), (16), (20)

Type of Test: Immunodot Comb

Manufacturer	Name of Test	Web Site	Useable Specimen	Perf. Ref.
Orgenics	HIV 1&2 BISPOT	http://www.orgenics.com	Serum, plasma	(16), (20)
Orgenics	IMMUNOCOMB II HIV 1&2 COMBFIRM	http://www.orgenics.com	Serum, plasma	
Orogenics	IMMUNOCOMB II HIV 1&2 SALIVA	http://www.orgenics.com	Saliva	
SPAN Diagnostics	Combaid -RS HIV1/HIV2	http://www.spanddiag.com	Whole blood	(5), (16),

				(17), (20)
Wierner Lab	DIA (Dot Immuno Assay) HIV 1+2	http://www.wierner-lab.com.ar	Serum, plasma	(5), (16)

HIV Performance References

1. Aidoo S et al. Suitability of a rapid immunochromatographic test for detection of antibodies to human immunodeficiency virus in Ghana, West Africa. *Journal of Clinical Microbiology* 39(7):2572-2575 (2001).
2. Amazu B et al. Evaluation of the rapid test "Determine HIV 1/2 (Trade Mark)" for the detection of HIV antibodies in Cameroon. *XIII International AIDS Conference Abstracts*.
<http://www.iac2000.org/abdetail.asp?ID=MoPeA2106> (8 April 2002).
3. Branson B et al. Performance of newer rapid tests for HIV antibody with whole blood and plasma. Centers for Disease Control and Prevention.
<http://www.medadvocates.org/conferences/Branson%20croi%20rt%20poster.pdf>.
4. Branson BM. Assessing diagnostic technologies marketed to less industrialized countries. *Journal of International Physicians in AIDS Care 2000*. February:28-30. Medical Advocates for Social Justice
http://www.medadvocates.org/journals/branson_paper.html.
5. Branson BM. HIV Rapid Testing: Challenges to Public Health. PowerPoint presentation.
<http://www.fda.gov/ohrms/dockets/ac/00/backgrd/3620b1p.pdf>.
6. Branson BM. Rapid test for HIV antibody. *AIDS Reviews 2000* 2:76-83. Medical Advocates for Social Justice
http://www.medadvocates.org/cdc/rapid_test.html (8 April 2002).
7. Burgess-Cassler A et al. A field test for the detection of antibodies to human immunodeficiency virus types 1 and 2 in serum or plasma. *Clinical and Diagnostic Laboratory Immunology* 3(4):480-482 (1996).

8. Kannangai R et al. Hospital-based evaluation of two rapid human immunodeficiency virus antibody screening tests. *Journal of Clinical Microbiology* 38(9):3445-3447 (2000).
9. Kassler WJ et al. Rapid HIV testing with same-day results: a field trial in Uganda. *International Journal of STD & AIDS* 9(3):134-138 (1998).
10. Kaur H et al. Evaluation of rapid kits for detection of HIV, HbsAg and HCV infection. *Indian Journal of Medical Sciences* 54(10):432-434 (2000).
11. Kline RL et al. Diagnosis and differentiation of HIV-1 and HIV-2 infection by two rapid assays in Nigeria. *Journal of Acquired Immune Deficiency Syndrome* 7(6):623-626 (1994).
12. Nookhai S et al. Evaluation of OraQuick HIV-1/2 rapid test and a rapid/simple test algorithm in a high-risk/high-prevalence population in Bangkok, Thailand. 8th Conference on Retroviruses and Opportunistic Infection.
<http://www.retroconference.org/2001/abstracts/abstracts/abstracts/232.htm>.
13. Palmer CJ et al. Field Evaluation of the Determine rapid human immunodeficiency virus diagnostic test in Honduras and the Dominican Republic. *Journal of Clinical Microbiology* 37(11):3698-3700 (1999).
14. Saville RD et al. Evaluation of two novel immunoassays designed to detect HIV antibodies in oral fluids. *Journal of Clinical Laboratory Analysis* 11(1):63-68 (1997).
15. Schramm W. A simple whole-blood test for detecting antibodies to human immunodeficiency virus. *Clinical and Diagnostic Laboratory Immunology* 5(2):263-265 (1998).
16. Simple and/or rapid HIV antibody test.
http://www.who.int/pht/blood_safety/hivtable2.html (5 April 2002).
17. Tribble DR et al. Comparative field evaluation of HIV rapid diagnostic assays using serum, urine, and oral mucosal transudate specimens. *Clinical and Diagnostic Virology* 7(3):127-132 (1997).

18. Urassa W et al. Evaluation of an alternative confirmatory strategy for the diagnosis of HIV infection in Dar Es Salaam, Tanzania, based on simple rapid assays. *Journal of Virological Methods* 100(1-2):115-120 (2002).
19. Walter Reed Army Institute of Research. "Technology for Rapid Transfusion Test Device (TTD) Testing. PowerPoint presentation. http://www.tricare.ods.mil/conferences/nato/downloads/4thplenarysession/2mi_chaelpresentation.ppt (April 2002).
20. World Health Organization. Operational Characteristics of Commercially Available Assays to Determine Antibodies to HIV-1 and/or HIV-2 in Human Sera, Report 11, Geneva, January 1999. WHO/BTS/99.1,UNAIDS/99.5.
21. Zaw M et al. Local evaluation of a rapid HIV assay for use in developing countries. *Tropical Medicine & International Health* 4(3):216-221 (1999).

See also:

[Diagnostic overview of rapid tests for HIV](#)
[Appropriate uses of rapid tests for HIV](#)

General Information

- [Technologies](#)
- [Appropriate Use](#)
- [Accuracy](#)

Tests for Specific Diseases

- [Hepatitis B](#)
- [HIV](#)
 - [Diagnostic overview](#)
 - [Appropriate use](#)
 - [Commercial availability](#)
- [Malaria](#)
- [Syphilis](#)